


SMARTACT *evo* SMARTACT PIN *evo*



The **evolution** goes on!
The unique features
of SMARTACT *evo* enable
the clinician to improve and
simplify the aspects regarding
fixing the membrane in the
precision of the surgical
operation, in access to all areas
of the oral cavity and in the
reduction of surgical times.



SMARTACT *evo* SMARTACT PIN *evo*




SMARTACT *evo*

The introduction of SMARTACT *evo* mechanical handpiece and the SMARTACT PIN *evo* titanium tacks revolutionizes the procedures and the traditional manual instruments used in regenerative oral surgery.

SMARTACT PIN *evo*

The new exclusive design of SMARTACT PIN *evo* allows absolute safety when inserting the tack on the handpiece. It gives great stability upon positioning on the bone and an axial insertion, facilitating the fixing of the membrane in every situation and hard-to-reach area. SMARTACT PIN *evo* SLIM is particularly suitable to use in dense mandibular bone.

The tacks are available in 4 different sizes, with retention and without retention, with flat head and standard head, to satisfy every clinical and aesthetic need. The tacks are sterile and ready to use.



courtesy of Dr. Luca De Stavola

Advantages

- Mechanical, ergonomic and balanced handpiece, small pedal
- Lightweight and easy to carry (max 980 gr.)
- Cable of three different lengths
- Absolute stability of the pin during positioning phase
- Silent with greater comfort for the patient
- Excellent access to every hard-to-reach area
- Easy to use and reduction of surgical times
- Different type of tacks according to bone density
- Tacks already sterile and ready to use
- Adjustable insertion force
- Reduction of risk factors
- Compressed air device

Ref. 5450-120 SMARTACT *evo* handpiece, pedal and cable 120 cm
Ref. 5450-150 SMARTACT *evo* handpiece, pedal and cable 150 cm
Ref. 5450-200 SMARTACT *evo* handpiece, pedal and cable 200 cm

 **Ref. 5060** SMARTACT PIN *evo* **STD red**
 pack of 3 tacks. Sterile

 **Ref. 5050** SMARTACT PIN *evo* **SLIM green**
 pack of 3 tacks. Sterile

 **Ref. 5100** SMARTACT PIN *evo* **FLAT blue**
 pack of 3 tacks. Sterile

 **Ref. 5090** SMARTACT PIN *evo* **SLIM FLAT white**
 pack of 3 tacks. Sterile




Dent-Med Materials
Specialists in Bone and Tissue Regeneration

T: 0226-360 150 info@dent-medmaterials.nl www.dent-medmaterials.nl

